The Swiss Instruments

- DE** Bedienungsanleitung
- EN** Operation Manual
- FR** Mode d'emploi
- IT** Istruzione d'uso
- ES** Instrucciones de uso

MD 11

Motor System for Implantology

DE

EN

FR

IT

ES

Vorwort

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Foreword

Congratulations on your purchase of a NOUVAG AG product. Thank you for the confidence shown in our products. Please consult the instruction manual for the use and maintenance of the device in order to ensure that it will function properly and efficiently for many years. You will find the conformity statement and list of authorized service representatives attached.

Federal Law restricts this device to sale by or on the order of a physician.

- **Please read instructions carefully before operating!**

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Table of contents

1	Product description	2
1.1	Intended use and operation	2
1.2	Contraindications	2
1.3	MD 11 technical data	2
1.4	Ambient conditions	2
1.5	Warranty coverage	2
2	Explanation of symbols	3
3	Safety information	3
3.1	EMC Manufacturer's Declaration of Conformity	3
3.2	Integrated peristaltic tube pump	3
3.3	Modification and misuse	3
3.4	Essential requirements	4
3.5	During use	4
4	Scope of delivery	4
5	Device overview	5
6	Startup	6
6.1	Device setup	6
6.2	Connection to the power supply	6
6.3	Device preparation	7
6.4	Assembly of external irrigation system	8
6.5	Assembly of internal and external irrigation system (optional)	9
7	7 Operation	10
7.1	Switching the device on and off	10
7.2	Overview: control elements on the operating panel	10
7.3	Overview: Standard display	11
7.4	Adjusting the programs	12
7.4.1	<i>Step 1: Selecting the handpiece or transmission ratio</i>	12
7.4.2	<i>Step 2: Calibrating handpieces or contra angles</i>	13
7.4.3	<i>Step 3: Setting the speed</i>	14
7.4.4	<i>Step 4: Setting the torque</i>	14
7.4.5	<i>Step 5: Setting the pump supply quantity</i>	14
7.5	Torque limit function, AL mode (Automatic Limiter)	15
7.6	Torque limit function, AS mode (Automatic Stopper)	15
7.7	Storing various programs	15
7.8	Configuration Menu	16
7.9	Operation using the Vario pedal	19
7.10	Functional check	20
7.10.1	<i>Electronic motor</i>	20
7.10.2	<i>Pump</i>	20
7.10.3	<i>Rotational direction of the electronic motor</i>	20
8	Cleaning, disinfection and sterilization	21
8.1	Control unit and pedal	21
8.2	Electronic motor 21	22
8.3	Tubing set, Art. No. 1706	23
8.4	Handpiece cradle	23
9	Maintenance	24
9.1	Replacing the control unit fuse	24
9.2	Safety inspections	24
10	Malfunctions and troubleshooting	25
11	Spare parts list with order numbers	28
12	Information on disposal	28

1 Product description

1.1 Intended use and operation

The MD 11 in combination with a motor and corresponding handpiece or contra angle (separate medical device) is used primarily in dental Implantology. The device is designed for drilling, milling and sawing bone as well as for screw insertion into bone. An integrated peristaltic pump is provided in order to cool the rotating instruments so that damage to tissue can be prevented.

The MD 11 may only be operated by qualified and trained personnel.

It is used specifically in Implantology for:

- Milling and drilling the implant bed
- Cutting the thread for the implant
- Screwing in the implant
- Removing the fixture mount
- Placing the cover screw

1.2 Contraindications

Not known.

1.3 MD 11 technical data

Voltage: ----- variable: 100 V~/ 115 V~/ 230 V~, 50–60 Hz
 Fuse, power supply: ----- 2 fuses T 3.15 AL 250 VAC
 Power consumption: ----- 120 VA
 Applied part: ----- Type B*
 Protection class: ----- Class I
 Dimensions (W x D x H): ----- 260 x 250 x 110 mm
 Net weight control unit: ----- 3.3 kg

Motor:

Motor coupling:----- Intra coupling ISO3964
 Motor speed: ----- 300 - 40,000 rpm
 Max. Motor torque: ----- 6 Ncm
 Motor weight: ----- 0.280 kg
 Motor cable length----- 2 m

Pedal:

IP code (pedal) ----- IPX8

Peristaltic pump:

Max. volumetric displacement of pump ----- 100 ml/min.

**Applied part of Type B is the instrument used with the MD 11.*


1.4 Ambient conditions

	Transport and storage:	Operation:
Relative humidity:	Max. 90%	Max. 80%
Temperature:	0°C to 60°C	10°C to 30°C
Atmospheric pressure:	700 hPa to 1060 hPa	800 hPa to 1060 hPa

1.5 Warranty coverage

Purchasing the MD 11 entitles you to a 1-year warranty. If you return the warranty card for registration within four weeks of the date of purchase, warranty coverage will be extended for a further **6 months**. Consumable parts are not covered by the warranty. Improper use or repair, or failure to observe these instructions, relieve us from any obligation arising from warranty provisions or other claims.

2 Explanation of symbols

	Important information		Autoclavable at 135°C for at least 5 min.
	Do not use if the packaging is damaged		Suitable for thermal disinfection
	Warning		Sterilized using ethylene oxide
	Manufacturer		Observe the instructions for use
1 min. on/ 3 min. off	The device is designed for intermittent duty operation at "1 min ON/3 min OFF" at 4 cycles, after 15 min. brake.		Electrical and electronic devices that have reached the end of their service life comprise hazardous waste and may not be disposed of together with household waste. Valid local disposal regulations apply.
	Type B applied part Applied parts is the instrument		Symbol indicating the serial number with the date of manufacture (year/month).
	Do not reuse		Symbol indicating the order number.
	Biohazard		Symbol indicating the lot number.
	Motor		Date of expiry
	Pedal		Protective ground
IPX8	Protection against contact and continual submerging.		Certified by the Canadian Standards Association (CSA)
	Date of manufacturing		CE symbol with notified body

EN

3 Safety information

Your safety, the safety of your team and of course that of your patients is very important to us. It is therefore essential to bear the following information in mind.

3.1 EMC Manufacturer's Declaration of Conformity

Please observe the information on electromagnetic compatibility provided with these operating instructions.

3.2 Integrated peristaltic tube pump

The integrated peristaltic pump is used to cool tissue in order to prevent damage to the tissue. It may only be operated with watery solutions, such as 0.9% Sodium Chloride irrigation solution (Art. No. 1707) or "Ringer" solution. Supplying medication using the integrated pump is expressly prohibited.

3.3 Modification and misuse

- Modification or manipulation of the MD 11 and its accessories is prohibited. The manufacturer is not liable for any damages resulting from unauthorized modifications or manipulations. The warranty will be canceled.
- Use of the MD 11 outside the indications described in Section 1.1 is prohibited. The user or operator is solely responsible for any such use.

3.4 Essential requirements

The MD 11 may only be operated by qualified and trained personnel.

Improper use or repair of the device and failure to observe these instructions relieve us from any obligation arising from warranty provisions or other claims.

The use of third-party products is the responsibility of the operator. Functionality and patient safety cannot be guaranteed with third-party accessories.

Prior to using the device, before startup, and before operation, the user must always ensure that the device and accessories are in good working order and are clean, sterile and operational.

Repairs may only be performed by authorized NOUVAG service technicians.

Use Nou-Clean spray for maintenance and care of the motor, handpieces and contra angles. Using other care products can result in malfunctions and/or cause the warranty to be revoked.

3.5 During use

The device is not sterile on delivery. All sterilizable parts must be sterilized before use (see Chapter 8 Cleaning, disinfection and sterilization).

At the choosing of the instrument the operator has to make sure it's biocompatible, according to EN ISO 10993.

Handpieces and contra angles may only be attached when the electronic motor is not running.

Never operate the clamping mechanism of the handpieces or contra angles while the system is running. This could result in instrument damage.

To prevent injury, never touch drill bits or burrs while they are still running.

Do not use device in the vicinity of flammable mixtures!

At operation at a patient, make sure to produce as less as possible friction heat. By excessive thermal influence on the tissue, the tissue easily gets necrotized. The heat generation is directly proportional to the speed delivery of the instrument and the contact pressure delivered by the instrument.

4 Scope of delivery

Art. No.	Description	Quantity
3335	MD 11 control unit-----	1 unit
1866nou	Vario pedal; IPX8; electronic -----	1 unit
2097	Electronic motor 21 incl. 2 m motor cable -----	1 unit
1706	Tube set, sterile, 2 m, single-use -----	1 unit
1873	Clip set (10 pieces) for tube set attachment at motor cable-----	1 unit
1881	Clip set (3 pieces) for tube set attachment at handpieces and contra angles-----	1 unit
1707	Irrigation fluid; 0.9% sodium chloride/water solution, 1 liter -----	1 unit
1770	Stand for irrigation fluid-----	1 unit
1170	Handpiece cradle-----	1 unit
1974	Spray nozzle for Nou-Clean spray; for care of electronic motor -----	1 unit
1958	Spray nozzle for Nou-Clean spray; for care of instruments -----	1 unit
31665	User Manual MD 11 on CD-ROM -----	1 unit


In line with regulations pertaining to hazardous materials, the following items are not delivered with the control unit and must be ordered separately:

1984	Nou-Clean maintenance and care spray -----	1 unit
------	--	--------

5 Device overview

Front view

- | | |
|---|---|
| 1. Pedal socket, device rear | 13. Tubing set |
| 2. Vario pedal | 14. Stand for irrigation fluid bottle |
| 3. Contra angle (optional) | 15. Drip chamber |
| 4. Clip for tubing set attachment at handpieces and contra angles | 16. Bleed valve |
| 5. Handpiece cradle | 17. Bottle with irrigation fluid |
| 6. Electronic motor | 18. Type plate with type designation, reference number, serial number, information on power supply and device fuse. |
| 7. Clip for tubing set attachment at motor cable | 19. Power plug socket |
| 8. Motor socket | 20. Main switch |
| 9. Operating panel | 21. Fuse compartment |
| 10. Display | 22. Spray nozzle for maintenance of electronic motor |
| 11. Release key for tubing set bracket | 23. Spray nozzle for maintenance of handpieces and contra angles. |
| 12. Peristaltic pump | |

Rear view

6 Startup

6.1 Device setup

- Place the MD 11 and all required accessories and instruments on an even, non-slip surface and make sure you have good access to all controls.
- Do not allow the operating range of the device (including the cable and contra angle) to be compromised by limiting factors.
- The system display must be fully visible at all times.
- The pedal must be placed within stepping distance between the patient and the surgeon.
- It must be explicitly ensured that no objects can fall on the pedal.
- The power plug at the rear of the device must be accessible at all times.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

6.2 Connection to the power supply

Before plugging the power cable into the power socket for the first time, you must check the supply voltage setting next to the power switch.

If the voltage shown does not correspond to the local mains voltage, the grey fuse holder must be set to the correct voltage:

- Unplug the power cable.
- Use a screwdriver to open the fuse slot.
- Remove the fuse holder.
- Remove the grey fuse holder and reinsert it so that the local mains voltage setting is shown in the small window.
- Slide the fuse holder back in and close the fuse slot.
- Check the mains voltage shown on the fuse slot.
- Plug the power cable back into the device.

In order to prevent the risk of electric shock, the device may only be connected to a power network with a PE protective ground conductor.

6.3 Device preparation

1. Sterilize the motor (the motor is not sterile on delivery). If the motor has already been sterilized: when removing the motor from the sterile packaging, ensure that the sterile packaging is not damaged and that the sterility indicator confirms sterility (if no sterility indicator is provided, the sterile packaging must at least show the date on which the shelf life of the sterile item is due to expire). 
2. Insert the stand for the irrigation fluid into the stand holder.
3. Plug the motor plug of the electronic motor into the motor socket.
4. Plug the pedal plug into the pedal socket at the rear of the control unit.
5. Attach the sterilized handpiece to the electronic motor. Press the handpiece firmly onto the electronic motor until it clicks into position and make sure it is secure by moving it slightly in the opposite direction.
6. Assemble the tubing set: Attach tubing set Art. No. 1706, for the cooling of the contra angle, like described in the following article.

Use only Nouvag tube set Art. No. 1706, otherwise the function cannot be guaranteed.

Check the expiry date of the tubing set and ensure that the packaging is not damaged. Using non-sterile tubing sets can result in serious infection and, in extreme cases, can be fatal.

When inserting the tubing set, watch the arrow on the cover of the pump compartment. It indicates the flow direction of the cooling liquid.

Do not regulate the amount of irrigation fluid using the roller clamp on the tubing set; with the MD 11, this is regulated instead using the integrated pump. For this reason, make sure to open the roller clamp as far as it will go (please refer to 7.4.5 Setting the pump supply quantity).

- A) Press release key for tubing set bracket (on top of the control unit) to open the pump.
- B) The compartment with the integrated tubing bracket opens.
- C) Place the tubing set into the tubing bracket provided in such a way that the part of the tubing set with the spike exits the pump towards the rear of the device. Make sure the tubing is secure.
- D) With the tubing set inserted, press the compartment downwards until it clicks into place.

7. Insert the spike at the end of the tubing set into the irrigation fluid bottle and hang the bottle onto the stand.
8. Open the roller clamp on the tubing set as far as it will go.
9. Open the bleed valve beneath the drip chamber.
10. Connect the control unit to the power socket.

Ensure that the operating voltage setting corresponds to the local mains voltage.


6.4 Assembly of external irrigation system

- A) Secure motor cable clips and handpiece clips to the irrigation tube.
- B) Connect tube set with the external cooling tube of the contra angle.
- C) Secure clips to the contra angle and to the motor cable.


If required, secure additional clips to the motor cable.

6.5 Assembly of internal and external irrigation system (optional)

- A) Graft a 6 cm tube piece (optional, Art. No.: 1773) on each of the branching tubes of the Y-Connector (optional, Art. No.: 1777).
- B) Connect one of the 6 cm tube pieces that branch out of the Y-Connector with the internal cooling tube (delivered with the contra angle, or optional, Art. No.: 1712).
- C) Secure small, grey clips and the bigger, white clips to the tube set.
- D) Secure bigger, white clips to the contra angle and the small, grey clips to the motor cable.
- E) Put the internal cooling tube with attached tube piece into the hole in the head of the contra angle and the drill bit (special drill bits are necessary with a hole for internal cooling).
- F) Push the spigot of the cooling tube into the hole in the shank of the clamping mechanism.
- G) Connect the 6 cm tube piece (optional, Art. No.: 1773) with the tube of the external cooling.

If needed, secure further clips to the cable of the motor or the contra angle.

EN


7 Operation

7.1 Switching the device on and off

The power switch "I/O" (at the rear) is used to switch the control unit on and off.

The device can be switched off at any time irrespective of any procedure for switching off the device.

7.2 Overview: control elements on the operating panel

A) **Display:** Shows the operating values (See chapter "7.3 Overview: Standard display").

B) **"Handpiece" key:** Selection of the 1:1 contra angle (Program 1).

- Bei pressing both keys at the same time "HP + PR", all Programs will be reset to factory default.

C) **"Program" key:** Selection of program 2 or 3 with the Contra Angle selected in the configuration menu. Select between 4 different Contra Angles, available in the configuration menu.

D) **"Speed" keys:**

- Restrict the maximum speed that can be selected using the pedal.
"▲": increases the maximum speed "▼": reduces the maximum speed
- By pressing both "Speed" keys at the same time "Speed ▲ + ▼" the handpiece calibration will be started (refer to "7.4.2 Calibrating handpieces").

E) **"Torque" keys:**

- Restricts the maximum torque.
"▲": increases the maximum torque "▼": reduces the maximum torque
- By pressing both "Torque" keys at the same time "Torque ▲ + ▼" the torque modes AL and AS are exchanged (see chapter "7.5 Torque limit function, AL mode (Automatic Limiter)" and chapter "7.6 Torque limit function, AS mode (Automatic Stopper)").

F) **"Pump" keys:**

- Changing the pump flow rate that can be supplied using the pedal.
"▲": increases the maximum supply quantity "▼": reduces the maximum supply quantity
- By pressing both "Pump" keys at the same time "Pump ▲ + ▼" the pump will be put on call, pressing again will switch it off.

7.3 Overview: Standard display

- A) **Information Line**
Information and error messages are displayed here. Display is illuminated red at error messages.
- B) **Program**
Shows selected number of program.
- C) **Pump**
The numerical value shows the pump flow rate in percent and the drop symbol together with the on/off indication shows if the pump is in stand-by mode or switched off.
- D) **Rotational direction of the motor**
The arrow indicates the rotational direction set for the motor. Change the rotation direction by pressing the button "FOR/REV" on the pedal.
- E) **Current speed**
At stand still of the device the maximum speed is displayed. If pedal is activated and motor begins to run the current speed is displayed.
- F) **Transmission ratio of contra angle**
Shows transmission ratio of used contra angle, e.g. 1:1, 16:1, 20:1, 32:1 or 70:1.
- G) **Speed range**
Shows speed range of the handpiece used.
- H) **Maximum torque**
Shows maximum torque setting.
- I) **Current torque**
Bar graph providing a graphical representation of the current torque. All bars active means max. torque reached.
- J) **AS/AL**
Indication of the Auto Stopper (AS) mode or the Automatic Limiter (AL) mode, (see chapter "7.5 Torque limit function, AL mode" and chapter "7.6 Torque limit function, AS mode").

The pump does not begin to operate until the motor has been activated by pressing the pedal.

7.4 Adjusting the programs

Values for operation settings depend on the connected handpiece or contra angle as well as the task to be performed.

7.4.1 Step 1: Selecting the handpiece or transmission ratio

Depending on the handpiece or contra angle attached to the motor, the corresponding transmission ratio has to be adapted accordingly by using the “HP” or the “PR” key.

Press the “**HP**” key to select Program 1 with the 1:1 contra angle
Press the “**PR**” key to select Program 2 or 3. The desired contra angle has to be selected in the configuration menu before you can use it.

Table of possible contra angles

Name of the handpieces/contra angles with transmission ratio	Display	Speed min. rpm	Speed max. rpm	Torque min. Ncm	Torque max. Ncm	AS-range (Factory def.) rpm	Limit AS-range rpm (*)
Drill handpiece, 1:1	1 : 1	300	40,000	1	6	No	–
Drill contra angle, 16:1	16 : 1	20	2400	5	27	up to 20	20 – 45 *
Drill contra angle, 20:1	20 : 1	15	1700	10	70	up to 20	15 – 45 *
Drill contra angle, 32:1	32 : 1	10	1000	10	55	up to 20	10 – 45 *
Drill contra angle, 70:1	70 : 1	5	600	10	55	up to 20	5 – 45 *

* The limitation of the AS-range (Automatic Stopper) can be adjusted in the Configuration Menu.

Handpieces or contra angles that don't belong to one's own assortment can be deactivated in the Configuration Menu. That makes finding the right handpiece in the handpiece list easier, because the list of available handpieces becomes shorter.

7.4.2 Step 2: Calibrating handpieces or contra angles

To make sure the displayed parameters correspond with the actual, measurable parameters of the handpiece or contra angle, it is recommended to calibrate each handpiece or contra angle on a regular basis.


It's a procedure as simple as it is important to guarantee safety and precision of each handpiece or contra angle being used.

After you take care of all prior preparations such as sterilization, maintenance and care of handpieces, device preparation and the selection of the handpiece of use, the calibration procedure can be performed.

The calibration of handpieces or contra angles guarantees accurate torque. Due to wear as well as varying lubrication of the handpieces or contra angles, and lack of maintenance and care, the distribution of torque can vary widely.

1. Select the handpiece, attached to the motor, by pressing the **"handpiece"** key. Make sure that this handpiece attached to the motor is also on the display.
2. Hold motor with mounted handpiece or contra angle in your hand at safe distance from your body.
3. Press both **"Speed"** keys at the same time (**Speed ▲ + Speed ▼**). Display indicates **"Calibrating handpiece XX"**.

4. Motor and handpiece start running and pass several speed cycles up to maximum speed.
5. After a tone emitted and a message displayed the calibration is finished. The display returns to normal mode.

If a handpiece is not working within the range of the calibration algorithms, even when it was cleaned and lubricated before, the display shows an error with a red backlit display, **„Handpiece XX is faulty“**. This indicates soiling, wear or a technical defect of the handpiece. These handpieces have to be repaired or replaced.

At the calibration procedure, the handpieces are tested on their torque delivery behavior. With the 1:1 Handpiece the control unit is additionally adjusted to altered conditions of the handpiece to keep the 1:1 Handpiece delivering the appropriate torque.

7.4.3 Step 3: Setting the speed

The possible speed range depends on the contra angle, which has to be activated in the Configuration Menu. The maximum speed within this speed range can be restricted to the required value. Using the pedal, the speed can be varied from the minimum speed up to the maximum speed as set.

Setting the speed:

Press the “Speed” keys “▲” to increase or “▼” to decrease the maximum speed. When key is pressed constantly the speeds will be shown in fast forward mode.

7.4.4 Step 4: Setting the torque

Once the speed has been selected, the torque can be set from the corresponding torque range. Depending on the speed, the torque mode AL or AS is applied.

For information on AL and AS mode, see Chapter “7.5 Speed reduction AL (Automatic Limiter)” and Chapter “7.6 Torque limit function: AS (Automatic Stopper)”


Press the “Torque” keys “▲” to increase or “▼” to decrease the maximum torque. By pressing the key continuously the torques are changed in fast forward mode.

The build up of torque is displayed as a bar graph. When maximum torque is reached all elements of the bar graph are full displayed.

7.4.5 Step 5: Setting the pump supply quantity

Press the “Pump” keys “▲” to increase or “▼” to decrease the pump supply quantity. By pressing the key continuously the values are fast forwarded.

Minimum and maximum of pump supply quantity as well as incremental steps can be adjusted in the Configuration Menu (refer to «7.8 Configuration Menu, Parameter level 2, Pump»).

To activate or deactivate the pump, press both “Pump” keys at the same time, “Pump ▲ + ▼”, or push the foot switch “PUMP” shortly.

50% OFF

7.5 Torque limit function, AL mode (Automatic Limiter)

The AL function limits the torque applied to the instrument, for example to prevent crack initiation and bone fracture.

The speed at the instrument remains constant until the preset torque is reached. If the applied force gets over this limit the speed will be reduced, if necessary down to a stop, but the torque remains constant. If the applied force is reduced the speed picks up again.

On the display this procedure is shown as a bar graph. The bars in the graph fill up to their full range as the torque increases to its preset value. When the bar graph shows full capacity the speed will be reduced. As soon as the force on the instrument is reduced the torque decreases and the speed starts to pick up again as shown on the display.

The AL mode is active at all speeds, except at the speeds where the AS mode takes effect.

7.6 Torque limit function, AS mode (Automatic Stopper)

The AS function limits the torque applied to the instrument. As soon as the preconfigured torque is reached, the electronic motor stops immediately. The electronic motor no longer generates any force. In order to restart the electronic motor, the treadplate must be released and pressed again.

On the display, the bar graph completes up to its full range, until the maximum torque is reached, then drops to zero.

The function **"AS"** is just active for certain contra angles and only in a certain speed range.

From the minimal speed up to the speeds shown at the right the AS-mode is activated automatically.
(refer to "7.8 Configuration Menu, Parameter level 2, AS-Zone")

Handpiece	16:1	20:1	32:1	70:1
Speed rpm:	20	20	20	20

The upper limit of the AS range can be adjusted in the configuration menu.

In the following speed ranges the AS-/AL-mode can also be switched manually by pressing both torque keys **"Torque ▲ + ▼"** at the same time.

Handpiece	16:1	20:1	32:1	70:1
From rpm:	20	15	10	5
To rpm:	45	45	45	45

7.7 Storing various programs

With the MD 11, up to 3 different settings can be set as fixed program (program 1 to program 3). Which program is currently active can be seen on the display.

At switch-off, the settings made by user are automatically saved. This includes the following parameters:

- Handpiece/Transmission ratio
- Maximal speed
- Maximal torque
- Pump On/Off
- Pump performance
- AS/AL mode

To change a program go to the specific parameter and change it. When the device is switched off all parameters are saved in that program.

When the MD 11 is switched on, the display shows the prior used program.

7.8 Configuration Menu

In the Configuration Menu the user can customize the device after his favor. The parameters are organized in different levels. The following information and parameters are stored in there.

- Software version
- Mainboard serial number
- Display illumination
- Operating hours MD 11
- Operating hours Motor
- Operating hours of irrigation pump
- Activation of available handpiece
- Speed limitation for each handpiece
- Range of action of AS mode
- Pump behavior
- Motor behavior
- Behavior of torque at reverse mode
- Reset to factory default

Be cautious when changing parameters. Unusual behavior of instruments while operating may provoke false reactions and jeopardize the patient. Every setting and the new behavior of the instrument has to be verified and tested.

Access to Configuration Menu:

- Press **"HP"** and **"Speed ▼"** key for about 3 sec. until you hear a tone. On the display the first position of the configuration menu shows up.

The arrows at the beginning of the lines indicate the operation in the Configuration Menu.

>Software
>Version V1.00

- Press **"HP"** or **"PR"** choose the desired parameter.
- Press **"Speed ▼"** to go to change mode.
- To change the value press **"HP"** for up or **"PR"** for down.
- To confirm the changes made, press **"Speed ▼"** for about one sec., until you hear a short tone.
- To abort the settings made press **"Speed ▼"** briefly, the setting returns to its previous value.

Exit from Configuration Menu:

- To leave the Configuration Menu press **"HP"** and **"PR"** for 3 sec., until you hear the long tone.

Parameter Level 1

Group/Parameter	Rights	Factory	Definition
Software/Version	read	VX.XXXX	Shows current software version
Hardware/Serial number MB	read	XXXXXXXXXX	Shows serial number of the main board
Backlight/brightness (0..10)	read/write	9	Display brightness, changeable: 0, ... , 10
Operating hours/Control Unit	read	0	Shows operating hours of Control Unit
Operating hours/Motor	read	0	Shows operating hours of motor 1
Operating hours/Pump	read	0	Shows operating hours of pump
Error memory/ 1 – 8	read	0	8 Error messages in chronological order.

Parameter Level 2

Values in Parameter Level 2 can only be adjusted after the password “9403” was entered. The password cannot be changed. Entering password: Press “Enter”, press “Select ▲ or ▼”, (for fast forward or backward keep key pressed).

EN

Handpiece activation	Name of handpiece on display	Choice	Factory default	Definition
Handpiece existing/HP 02	16 : 1	yes/no	no	Only one handpiece or contra angle can be selected.
Handpiece existing/HP 03	20 : 1	yes/no	yes	
Handpiece existing/HP 04	32 : 1	yes/no	no	
Handpiece existing/HP 05	70 : 1	yes/no	no	

Handpiece Max. Speed	Name of handpiece on the display	Speed range rpm	Factory default	Definition
Handpiece max speed/HP 01	1 : 1	300 – 40,000	40,000	Limit the maximum speed of your handpieces according to your own experience.
Handpiece max speed/HP 02	16 : 1	20 – 2400	2400	
Handpiece max speed/HP 03	20 : 1	15 – 1700	1700	
Handpiece max speed/HP 04	32 : 1	10 – 1000	1000	
Handpiece max speed/HP 05	70 : 1	5 – 600	600	

AS-Zone for handpieces	Handpiece on display	Speedrange ... to ... rpm	Factory default	Definition
Handpiece AS-Mode/HP 06	16 : 1	20 – 45	20	Effective AS-range 20 – 45 rpm
Handpiece AS-Mode/HP 07	20 : 1	15 – 45	20	Effective AS-range 15 – 45 rpm
Handpiece AS-Mode/HP 08	32 : 1	10 – 45	20	Effective AS-range 10 – 45 rpm
Handpiece AS-Mode/HP 09	70 : 1	5 – 45	20	Effective AS-range 5 – 45 rpm

Pump parameters	Range	Factory default	Definition
Pump/Backwards turn mode variable	No/Yes	Yes	The pressures in the tube set vary according to the pump speed. In “variable mode” this is considered, to prevent of spilling when pump is switched off in backwards mode.
Pump/Way backwards	1 – 100%	25%	Specify how far the pump turns backwards
Pump/Speed backwards	10 – 50%	33%	Specify how fast the pump has to turn backwards to prevent of spilling after the switch off of the motor.
Pump/Range 1 increment	1 – 10%	5%	Adjustment steps in section 1.
Pump/Range 1 end	5 – 50%	10%	Set the range where section 1 is active
Pump/Range 2 increment	1 – 10%	5%	Adjustment steps in section 2
Pump/Range 2 end	10 – 90%	50%	Set the range where section 2 is active.
Pump/Range 3 increment	1 – 10%	10%	Adjustment steps in section 3
Pump/Range 3 end	20 – 100%	100%	Set the range where section 3 is active.

Motor type 21, 40,000 rpm	Range	Factory default	Definition
Motor type 21, 40,000 rpm/Min. speed	300 – 5000 rpm	300 rpm	Set the min. speed of the motor.
Motor type 21, 40,000 rpm/Max. speed	5000 – 40,000 rpm	40,000 rpm	Set max. speed of the motor.
Motor type 21, 40,000 rpm /Start ramp	1 – 1000 ms/10,000 rpm	100 ms	Set acceleration time to 10'000 rpm
Motor type 21, 40,000 rpm /Stop ramp	1 – 1000 ms/10,000 rpm	50 ms	Set breaking time from 10'000 – 0 rpm

The maximum torque of the motor is boosted at the starting procedure of the reverse rotation to disengage stuck implants or screws.

Reverse torque	Range	Factory default	Definition
Reverse torque/Increase	5 – 30%	25%	Increase of the selected torque in reverse rotation.
Reverse torque/Increase time	100 – 2000 ms	500 ms	Time during which torque is raised.

Resetting to factory default	Choice	Factory default	Definition
Default value/Set default value	Yes/No	–	Resetting all parameters of the Configuration Menu to factory default.


- **Attention:** Upon resetting to factory default all parameters will appear with factory default values (except date and time).

Exit from Configuration Menu:

- To leave the Configuration Menu press **“HP”** and **“PR”** for 3 sec., until you hear the long tone.

7.9 Operation using the Vario pedal

EN

1. **Key "PUMP":**
Pressing the key briefly: switches the pump on or off (see information on display).
Pressing the key longer: increases the pump speed (see information on display).
2. **Key "FOR/REV":**
Pressing the key briefly: switches the direction of revolution. (see information on display).
3. **Treadplate**
With the treadplate pedal the motor speed is variably adjusted and the pump is activated.

Treadplate...	Motor:	Pump:
... not pressed	Motor off	Pump off
... pressed gently	Motor runs slowly	Pump on, if pump "On" displayed (speed as set on the control unit)
... pressed all the way down	Motor runs at maximum speed (speed as set on the control unit)	Pump on, if pump "On" displayed (speed as set on the control unit)

For safety reasons, the unit can only be operated by pedal.

7.10 Functional check

Prior to MD 11 startup or use of accessory equipment, the user must always ensure that each individual component is in good working order, free from defects, and is clean, sterile and operational. All inscriptions on the device and its accessories must be readable and there must be no loose parts in the device. Once the device is switched on, the most recent settings entered appear on the display.

7.10.1 Electronic motor

Use the **“Speed”** selection keys to set the speed of the electronic motor to its maximum. With the 1:1 Contra Angle attached, it is 40,000 rpm. With all the other Contra Angles the maximum Speed is lower, according to the transmission ratio of the Contra Angel. Press the pedal treadplate; the electronic motor starts up and accelerates to up to 40,000 rpm. When the treadplate is released, the electronic motor slows down again.

- The electronic motor is designed for an intermittent duty operation at maximum speed of “1 min ON/3 min OFF” at 4 cycles. Brake after 15 min. Otherwise the electronic motor can be damaged by excessive heat generation, and by touching the motors hull it could cause serious burnings.
- The motor ventilation slots must be kept clear in order to prevent the motor temperature from becoming excessive.

7.10.2 Pump

Press the key **“PUMP”** on the pedal briefly; the peristaltic pump is switched on, which is shown on the display by the symbol of a drop. Press the pedal treadplate; the peristaltic pump and the electronic motor start up. Water sprays from the irrigation needle on the contra angle.

7.10.3 Rotational direction of the electronic motor

Press the key **“FOR/REV”** on the pedal briefly; the rotational direction of the electronic motor changes. Press the pedal treadplate; the electronic motor rotates to the left and a continuous tone is emitted. Release the treadplate; the electronic motor ceases to operate and the tone is no longer heard. By pressing the motor key again, the rotational direction is switched back to right rotation, which is shown on the display by the symbol of a changing arrow.

8 Cleaning, disinfection and sterilization

The following points in particular are important with regard to caring for the material:

- Perform cleaning, disinfection and sterilization after every treatment.
- Always autoclave the material in sterilization packaging.
- Make sure that sterilization packaging is no more than 80 % full.

- Always autoclave the material at 135°C for at least 5 minutes.
- If sterilized material is not used immediately, the material packaging must be labeled with the sterilization date.
- Nouvag AG recommends including a sterility indicator.

EN



8.1 Control unit and pedal

Wipe the outside using micro-biologically tested surface disinfectant or a 70% isopropyl alcohol solution. The front plate of the control unit is sealed for this purpose and can be wiped clean.

8.2 Electronic motor 21

Reprocessing restrictions	Frequent reprocessing has only a limited impact on the electronic motor. The end of the product service life is normally determined by wear and damage through use. Electronic motor 21 is designed for 250 sterilization cycles.
----------------------------------	---

INSTRUCTIONS	
At the location of use	Remove surface soiling with a disposable cloth/paper towel.
Storage and transport	No special requirements. Due to the risk of drying and corrosion, reprocessing must be performed without undue delay.
Preparation for cleaning	<p>Remove soiling from the electronic motor with a disposable cloth/paper towel. Unscrew the motor cap, remove the cable and unscrew the handpiece carrier.</p> 
Automatic cleaning and disinfection	<p>Equipment: Washer-disinfector with a special load carrier that ensures motors are connected to the washer-disinfector and channels are rinsed. Rinse the motor from the front. Only use a neutral cleaning agent for this purpose.</p> <p>Place the electronic motor in the load carrier (ensure that the channels can be rinsed). Place the motor cap and cable, and handpiece holder, in the basket.</p> <p>Set a cleaning cycle that offers sufficient cleaning and rinsing. Perform the final rinse with fully deionized water.</p> <p>Perform a 10-minute rinse cycle at 93°C to facilitate thermal disinfection.</p> <p>When removing the electronic motor, check the motor cap and cable, and handpiece holder, to verify whether soiling is still visible in the gaps and the grooves. If necessary, repeat the cycle or clean manually.</p>
Manual cleaning	<p>Equipment: Neutral cleaning agent, soft brush, running demineralized water (< 38°C)</p> <p>Procedure:</p> <p>Rinse off and brush away surface soiling on the electronic motor, motor cap and cable, and handpiece holder.</p> <p>Use a brush to apply cleaning agent to all surfaces and gaps.</p> <p>Rinse the electronic motor, motor cap and cable, and handpiece holder thoroughly under running water.</p>
Manual disinfection	For manual disinfection, wipe the electronic motor, the motor cap and the plug and cable underneath it, and the handpiece holder, with a certified non-chlorine disinfectant.
Drying	If a drying program is not provided by the washer-disinfector, the electronic motor must be dried manually or in a drying cabinet. The handpiece holder must then be screwed back onto the motor.

Inspection and maintenance  <p>Art. No. 1974</p>	<p>Perform a visual inspection to check for damage, corrosion and wear. Connect Motor with the spray nozzle (Art. No. 1974) by screwing it on and spray motor for about 3 seconds and briefly wipe with a moist cloth (see instructions on spray can). Once thorough spraying has been completed, screw the cable and motor cap back onto the electronic motor.</p> 
Packaging	<p>Individual: Pack the electronic motor in individual packaging for sterile items. The bag must be large enough to ensure that the seal is not subject to strain. Nouvag AG recommends including a sterility indicator.</p> <p>Sets: Sort electronic motors on trays intended for this purpose or place them on all-purpose sterilization trays.</p>
Sterilization 	<p>Autoclave in the vacuum autoclave (Class B or S as per EN 13060) at 135°C for at least 5 minutes*. When sterilizing several instruments during one sterilization cycle, do not exceed the maximum sterilizer load. A drying cycle must be added in the case of autoclaves without a post-vacuum function. Allow the electronic motor to dry in the bag for at least 1 hour at room temperature with the paper side facing upwards.</p> <p>* Temperature exposure times are based on country-specific guidelines and standards. Maximum exposure time: 25 minutes.</p>
Storage	<p>If the sterilized electronic motor is not used immediately after sterilization, the material packaging must be labeled with the sterilization date. Including a sterility indicator is recommended.</p>

The effectiveness of the sterilization instructions provided above for electronic motor 21 and the handpiece has been validated by Nouvag AG. The user is responsible for ensuring that the sterilization procedure performed achieves the required result. This requires validation and routine monitoring of the procedure. The staff member who completes the procedure bears sole responsibility for any deviation on his part from the instructions provided. Deviations necessitate revalidation of the effectiveness of the procedure as well as of the technical resilience of the motor with regard to the modified sterilization process.

8.3 Tubing set, Art. No. 1706

- Single-use tubing set 1706 may not be reused.
- Tubing sets must be disposed of properly after use!
- Do not use tubing sets when packaging is already open or damaged!
- Do not use tubing sets when expiry date has passed.
- Use only Nouvag tubing sets with Art. No. 1706!

Sterility cannot be guaranteed by reusing and re-sterilization of tubing sets. The characteristics of the material change in a manner that can result in failure of the system. This may result in serious infections or even patient death worst-case.

8.4 Handpiece cradle


Soiled handpiece cradles are cleaned using a neutral cleaning agent and then sterilized in accordance with the same instructions as for electronic motor 21.

9 Maintenance

9.1 Replacing the control unit fuse

Users can replace faulty control unit fuses themselves. These are located at the rear of the device in the fuse slot beside the power switch:

- Unplug the power plug.
- Open the fuse slot using a screw driver.
- Replace the faulty fuse T 3.15 AL 250 V AC.
- Slide the fuse holder back in and close the fuse slot.
- Check the mains voltage shown on the fuse slot.
- Plug in the power plug again.

1. Fuse slot locking mechanism
2. Display window for voltage setting
3. Fuse slot
4. Fuse 1
5. Fuse 2

9.2 Safety inspections

The performance of safety inspections on medical devices is required by law in several countries. The safety inspection is a regular safety check that is compulsory for those operating medical devices. The objective of this measure is to ensure that device defects and risks to patients, users or third parties are identified in time.

An inspection interval of **2 years** applies to the MD 11.

NOUVAG AG offers a safety inspection service for its customers. Addresses can be found in the appendix of this operation manual under "Service centers". For further information please contact our technical service department.

Further international service centers are listed on the Nouvag website:

www.nouvag.com > Service > Service centers

10 Malfunctions and troubleshooting

Malfunction	Cause	Solution	Refer to operating instructions
Device is not operational	Control unit is not switched on	Set the power switch "I/O" to "I"	7.1 Switching the device on and off
	Power connection not established	Connect the control unit to the mains power supply	6.2 Connection to the power supply
	Incorrect operating voltage	Check the mains voltage	6.2 Connection to the power supply
	Faulty fuse	Replace the fuse	9.1 Replacing the control unit fuse
Motor does not run	Motor not switched on	Switch on the motor using the treadplate	7.9 Operation using the Vario pedal
	Motor not connected	Connect the motor cable to the control unit	5.0 Device overview 6.2 Connection to the power supply
	Handpiece or contra angle not correctly assembled	Press the handpiece firmly onto the electronic motor until it clicks into position and check that it is secure by moving it slightly in the opposite direction.	6.3 Device preparation
No irrigation fluid for instrument	Peristaltic pump not switched on	Switch on the peristaltic pump	7.9 Operation using the Vario pedal
	Tubing set incorrectly inserted	Insert tubing set correctly (note the direction)	6.3 Device preparation
	Tubing set clogged/crusted matter visible	Replace the tubing set	6.3 Device preparation
	Bottle with sodium chloride solution not ventilated	Open the ventilation filter at the drip chamber	6.3 Device preparation
	Tubing set is dripping	Replace the tubing set	6.3 Device preparation
	Roller clamp of tubing set is closed	Open roller clamp all the way	6.3 Device preparation
Pedal is not operational	Pedal not connected	Connect the pedal to the control unit	6.3 Device preparation
	Incorrect operation	Check operating instructions	7.9 Operation using the Vario pedal

If a fault cannot be rectified, please contact your supplier or an authorized service centre. The addresses are provided on the last page of the operating instructions.

EN

MD11, Error-messages on Display

Error-messages/Error Code	Cause	Solution
Basic Initialization/ W00	First Initialization	
Set default value/ W01	Parameter reset to default value	
Memory error/ E02	System Error	Send Control Unit to Service Center.
Handling error/ E03	System Error	Send Control Unit to Service Center.
Program SW error/ E04	System Error	Send Control Unit to Service Center.
User Config SW error/ E05	System Error	Send Control Unit to Service Center.
Display error/ E06	System Error	Send Control Unit to Service Center.
Pump error/ E07	System Error	Send Control Unit to Service Center.
Storing factory settings/ User Config & Program	Message while default values of parameters and programs saved to the NOU-dongle.	
Storing factory settings/ Program	Message while default values of programs are stored.	
Pedal not connected/ E10	a) Pedal not plugged in. b) Plug or cable is defective.	a) Plug in Pedal b) Send Control Unit and pedal to Service Center.
Pedal test mode/ W11	Pedal test mode switched on	Switch off device for 5 seconds, than switch on again.
Keyboard test mode/ W12	Keyboard test mode switched on	Switch off device for 5 seconds, than switch on again.
No motor connected/ E13	a) No motor connected b) Motor, motor cable, motor plug or Control Unit are is defective	a) Plug in motor b) Send motor and Control Unit to Service Center.
Unknown motor/ E16	a) Wrong motor plugged in. b) Right motor connected, but motor, motor cable, motor plug or Control Unit is defective.	a) Plug in correct motor b) Send motor and Control Unit to Service Center.
Pump is open/ E20	Motor is not working when pump compartment is open to prevent of injuries.	Close pump compartment.
Motor or pump test mode/ W21	Motor or pump test mode is switched on.	Switch off device for 5 seconds, than switch on again.
AS-mode torque reached	When max. torque setting is reached in AS-mode, this message is displayed.	Release pedal and restart motor by pressing pedal again.
Pedal locked/ W26, pedal release	If pedal is pressed at switch on procedure, pedal will not work.	Release pedal for one second.
Handpiece XX is faulty/ E29	The handpiece was exposed to excessive torque, while calibrating or testing.	<ul style="list-style-type: none"> - Clean handpiece and spray it thoroughly with Nou-Clean spray. - If message still displayed after testing, send handpiece to Service Center.

Handpiece XX is Ok!	Tested handpiece is OK.	
Calibrating HPXX is Ok!	Calibrated handpiece is OK.	
Testing the handpiece XX	Handpiece is testing.	
Calibrating handpiece XX!	Handpiece is calibrating.	
Nou-dongle is plugged in	Message is displayed for 1 second after NOU-dongle was plugged in.	

11 Spare parts list with order numbers

Accessories	Art. no.
Clip set large CL, for attachment to the handpiece, package with 3 pieces	1881
Clip set motor cable, for attachment to the motor cable, package with 10 pieces	1873
Single-use tubing set, 2 m, sterile, 10-unit pack	1706
Y-connector, for branching tube set for internal and external cooling	1777
Internal cooling nozzle for drilling handpieces with latch system	1712
Internal cooling nozzle clip for drilling handpieces with push button system	39116
Irrigation fluid, 0.9% sodium chloride, 1 liter	1707
Nou-Clean clean and care spray	1984
Spray nozzle attachment for surgical instrument e-coupling	1958
Spray nozzle attachment for electronic motor 21	1974
Surgery Handpiece 1:1, length 90 mm, with attachable spray nozzle	1710
Surgery Handpiece 1:1, length 105 mm, with attachable spray nozzle	1950
Surgery Contra Angle 1:1, length 125 mm, with attachable spray nozzle	1960
Drilling Contra Angle 1:1, with INTRA EN3964 coupling, max. torque 10 Ncm	5050
Drilling Contra Angle 16:1, with INTRA EN3964 coupling, max. torque 30 Ncm	5060
Drilling Contra Angle 20:1, with INTRA EN3964 coupling, max. torque 80 Ncm	5053
Drilling Contra Angle 20:1 LED, with INTRA EN3964 coupling, max. torque 80 Ncm	5052
Drilling Contra Angle 32:1, with INTRA EN3964 coupling, max. torque 55 Ncm	5054
Drilling Contra Angle 70:1, with INTRA EN3964 coupling, max. torque 55 Ncm	5065

To order additional parts, please contact our customer service department.

Manual MD11----- 31665

At Nuvag AG the manual is delivered as a PDF-File on CD-ROM, together with the according device. If you prefer a printed Version of it we will send you one, free of charge, by postal delivery. If you cannot find the manual anymore, we'll send you a spare one in PDF-Format by E-Mail.

12 Information on disposal

When disposing of the device, device parts and accessories, the regulations prescribed by law must be observed.

Do not dispose of devices with household waste!

To ensure environmental protection, old devices can be returned to the dealer or manufacturer.

Motors that have reached the end of their service life may not be disposed of with household waste. Motors must be sterilized before disposal. Please observe currently valid national disposal regulations for infectious waste.

Contaminated single-use tubing sets are subject to specific disposal requirements. Please observe currently valid national disposal regulations for infectious waste.

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Electromagnetic compatibility (EMC)

Remark:

The **Product** subsequently referred to herein always denotes the MD 11.

Changes or modifications to this product not expressly approved by the manufacturer may result in increased emissions or decreased immunity performance of the product and could cause EMC issues with this or other equipment. This product is designed and tested to comply with applicable regulations regarding EMC and shall be installed and put into service according to the EMC information stated as follows.

WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the product may cause unexpected or adverse operation.

WARNING

The product shall not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the product shall be tested to verify normal operation in the configuration in which it is being used.

Compliant Cables and Accessories

WARNING

The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the product.

The table below lists cables, transducers, and other applicable accessories for which the manufacturer claims EMC compliance.

NOTE: Any supplied accessories that do not affect EMC compliance are not listed.

Description	Length max.
Electronic motor 21	2.9m
Vario Footpedal IPX8	2.9m

Guidance and manufacturer's declaration – electromagnetic emissions

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	< 5 % U_T (> 95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (> 95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Product requires continued operation during power mains interruptions, it is recommended that the Product be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity for not life support equipment			
The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.			
Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
			<p>Portable and mobile RF communications equipment should be used no closer to any part of the Product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p>
Conducted RF IEC 61000-4-6	3 V rms 150 kHz to 80 MHz outside ISM bands	10 V rms 150 kHz to 80 MHz outside ISM bands	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz
			<p>where P is the maximum output power rating in the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a	Fixed strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Product is used exceeds the applicable RF compliance level above, the Product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Product.		
b	over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended separation distances between portable and mobile RF communications equipment and the not life support equipment			
The Product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Product as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the higher frequency range applies.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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For all Nouvag servicecenters worldwide please check:

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Nouvag AG • St.Gallerstr. 23-25 • CH-9403 Goldach
Nouvag GmbH • Schulthaißstrasse 15 • D-78462 Konstanz
Nouvag USA Inc. • 18058 Albyn Court • Lake Hughes, CA 93532 • USA